

MANUKA WOUND CARE



HONEY ENHANCES WOUND CARE FOR VETERINARIANS.

CLINICALLY PROVEN. ALL NATURAL. EASY TO USE.


JorVet is proud to offer our 100% active medical-grade Manuka Honey Wound Care Dressings that help maintain a moist wound environment conducive to enhancing wound healing. The high sugar levels of our medical-grade Manuka honey result in osmotic pressure which promotes autolytic debridement of necrotic tissue and rapid reduction in wound malodor.

Extensive research has determined that Manuka honey contains synergistic factors that yield properties to maintain a balanced, moist wound environment. Our Manuka honey contains high levels of Methylglyoxal (MGO) and is standardized for anti-bacterial efficacy and sterilized for safety.

Under the supervision of a veterinary professional, Manuka Honey Wound Care Dressings may be used for:

- Trauma Wounds
- Surgical Wounds
- Ulcerated Wounds
- Infected Wounds
- 1st & 2nd Degree Burns
- Lacerations
- Abrasions
- Hard to Heal Wounds

THE LOW pH OF MANUKA HONEY CREATES AN ACIDIC ENVIRONMENT CONDUCTIVE TO IMPROVING WOUND HEALING.



Manuka Honey Dressings to match all types of wound needs.


MANUKA FILL®

- Honey gel in a tube
- All-purpose wound dressing
- Perfect for deep, tunneling and hard to dress areas



1.5 oz

0.5 oz




MANUKA IG®



4" x 5"
10cm x 12.5cm

- Non-adherent single layer gauze sheet with Manuka Honey
- Light absorbency
- Non-tacky surface




MANUKA IG® MAX



4" x 5"
10cm x 12.5cm

- Non-adherent 3 layer gauze sheet with Manuka Honey
- Light to moderate absorbency
- Easy to apply and remove from wound
- Can be trimmed to smaller sections and remainder stored safely in packaging




MANUKA FOAM® AIR



4" x 4"
10cm x 10cm

- Non-adherent absorbent foam pad with Manuka Honey
- Moderate to high exudate absorbency
- Easy to apply and remove from wound



All Manuka Honey Wound Care Dressings protect the wound while creating a moisture balance around the wound that is ideal for healing.

VETERINARY MANUKA WOUND CARE LINE

Description	Product Code	Pkg/ Unit
MANUKA FILL® Sterile Manuka Honey Wound Dressing		
15g Tube	J1253	10
42.5g Tube	J1253A	10
MANUKA IG® Manuka Honey Impregnated Gauze		
4" (10cm) x 5" (12.5cm)	J1251A	10
MANUKA IG® MAX Manuka Honey Impregnated Customized Gauze		
4" (10cm) x 5" (12.5cm)	J1252A	10
MANUKA FOAM® AIR Manuka Honeycomb Foam Dressing		
4" (10cm) x 4" (10cm)	J1254	10



Jorgensen Laboratories
1450 N. Van Buren Ave.
Loveland, CO 80538
(800) 525-5614
www.jorvet.com
info@jorvet.com